

CANIS 5A PORTABLE X-RAY MACHINE



FEATURES

- ❏ Compact and Lightweight
- ❏ Easy to Operate, Assemble and Transport
- ❏ Powerful 5.3 kW Generator
- ❏ Collimator Built-in Laser Ranging Light Source
- ❏ Super Affordable



For More Information Contact
1-305-463-9447

Biomedical Equipment Corp. 10421 NW 28TH ST STE D102,
DORAL, FL 33172 | E: info@4biomed.com | W: www.4biomed.com

SPECIFICATIONS

Model

Canis 5A

System Requirements

Electric Source (V)

AC 210 - 230 / AC 100 - 120

Power (KW)

5

Generator

kV Range (KVp)

40-125

mAs Range (mAs)

0.1-100

Ranging Mode

Ultrasonic

Communication Mode

RS232/Bluetooth

Exposure Mode

Wired/Wireless

X-ray Tube

Focal Spots (mm)

0.6/1.8

Anode Angle (°)

15


Mechanics Parameters

Size (mm³)

370 x 350 x 230

Weight (kg)

16



MOUNTING OPTIONS



Biomedical Equipment Corp. 10421 NW 28TH ST STE D102,
DORAL, FL 33172 | E: info@4biomed.com | W: www.4biomed.com